HemoCue® Plasma/Low Hb

Operating manual

Manual de funcionamiento

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EN

Table of Contents

HemoCue® Plasma/Low Hb System.....	3
Components.....	4
Start-up	6
Measuring Sample material	8
Maintenance.....	12
Troubleshooting Guide.....	14
Specifications.....	23

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HemoCue® Plasma/Low Hb System

EN

Thank you for choosing the HemoCue Plasma/Low Hb System. The system is intended for quantitative determination of low levels of hemoglobin in plasma and serum specimens, aqueous solutions, or supernatant from erythrocyte suspensions using a specially designed photometer, the HemoCue® Plasma/Low Hb Photometer and specially designed microcuvettes, the HemoCue® Plasma/Low Hb Microcuvettes.

The HemoCue Plasma/Low Hb Photometer is only to be used with HemoCue Plasma/Low Hb Microcuvettes. The HemoCue Plasma/Low Hb Microcuvettes are for *in vitro* diagnostic use only.

Read and follow this operating manual and the HemoCue Plasma/Low Hb Microcuvettes package insert to attain optimum performance and safety. Any other use of the system than recommended by the manufacturer may impair the safety.

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
Materials Provided

1. HemoCue Plasma/Low Hb Photometer*
2. Power adapter**
3. HemoCue Plasma/Low Hb Operating Manual

Open the carton and lift out the photometer and other components.
If no AC power is available, use 5 alkaline type AA batteries. On the bottom of the photometer there is a lid covering the battery compartment. Press the flap to remove the lid. Place the batteries in the battery compartment and replace the lid.
Consult local environmental authorities for proper disposal of batteries.

* Do not open the cover of the photometer. The warranty is void if the cover of the photometer has been opened.

** ⚠ Only use the power adapter listed under *Specifications*.

For information about HemoCue Plasma/Low Hb Microcuvettes please contact the HemoCue distributor.

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



Start-up
Activación
起動

Make sure the photometer is placed horizontally and on a stable surface.

 Only use the power adapter approved by HemoCue, as listed under *Specifications*.

1. If AC power is available, use the power adapter.
2. If no power is available, insert 5 alkaline type AA batteries.
3. Turn the switch to the 'ON' position (3a). All symbols appear on the display (3b), and after approximately 10 seconds the letter 'LHb' will be displayed (3c).
4. Pull out the cuvette holder to its loading position. This point, which should not be exceeded, is easily established by paying attention to a distinct stop.
5. After approximately 15 seconds the indication 'READY' appears on the display together with three flashing dashes. The photometer is ready for measurement.

Turn the photometer off by switching power switch to the 'OFF' position.

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



④

⑤

- ⚠ Always wear protective gloves. Handle blood with care, as it may be infectious. Follow local safety procedures for disposal of used microcuvettes.

Visually turbid samples should be filtered before analysis using a 0.2 µm filter.

1. To perform a test, the cuvette holder should be in its loading position. The display will show three flashing dashes and 'READY'. Remove a HemoCue Plasma/ Low Hb Microcuvette from the vial.
2. Place a drop of well mixed sample onto a hydrophobic surface using a suitable transfer device.
3. Fill the microcuvette in one continuous process. **Do not refill.**
4. Wipe off excess sample from the outside of the microcuvette with a clean, lint-free wipe. Do not touch the filling end of the microcuvette. If a second sample is to be taken, fill a new microcuvette from a new drop of sample. This should not be done until the measurement of the first sample is completed.
5. Visually inspect the microcuvette. If the microcuvette is not completely filled with sample or there are air bubbles, discard and fill a new microcuvette. Small bubbles around the edge can be ignored.

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




The cuvette holder should be cleaned after each day of use. Make sure that the photometer is turned off and that the display is blank.

1. Pull the cuvette holder completely out. Clean the cuvette holder with alcohol (20–70%) or mild detergent. The cover may be cleaned with alcohol (20–70%) or a mild detergent.
2. If the optical parts become dirty, an error code will be displayed. Before cleaning the optical parts, the cuvette holder should be completely removed as described above in step 1. To clean the optical parts inside the photometer, push a HemoCue Cleaner swab into the opening of the optic unit. For optimal cleaning push the cleaner only as far as until the white part of the cleaner swab is inside the optic unit. Move the cleaner back and forth 5–10 times. If the cleaner swab becomes stained, repeat with a new cleaner. Wait 15 minutes before replacing the cuvette holder.

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Troubleshooting Guide

If you are unable to resolve the problem by following this Troubleshooting Guide, please contact the local HemoCue distributor or HemoCue AB. Prior to service or disposal, the photometer should be cleaned as recommended under *Service and disposal* section. The photometer has no serviceable parts. Do not open the cover of the photometer. The warranty is void if the cover of the photometer has been opened.

Symptom	Explanation	Action
The photometer shows an error code	May be a temporary fault.	Turn off the photometer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
ERROR 900	No stable endpoint of the measurement is found within the time range. 1. The microcuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiration date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The photometer needs service. Contact the distributor.
ERROR 903	1. Disturbances on main power supply. 2. The optronic unit is out of order.	1. Change the wall socket or use battery power. 2. The photometer needs service. Contact the distributor.
ERROR 901–905	1. Dirty optical parts or faulty electronics or optronic unit.	1a. Turn off the photometer and clean the optical parts as described in the <i>Maintenance</i> section. 1b. The photometer needs service. Contact the distributor.
ERROR 906	1. Unstable blank value. The photometer might be cold.	1. Turn off the photometer and allow it to reach operating temperature. If the problem continues, the photometer needs service. Contact the distributor.
ERROR 907	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the photometer and replace the batteries, 5 alkaline type AA. 1b. Use the power adapter.
ERROR 908	The absorbance is too high. 1. Light blocking item in the cuvette holder.	1a. Check that the photometer and microcuvettes are used according to the HemoCue Plasma/Low Hb operating manual and package insert. 1b. The photometer needs service. Contact the distributor.
ERROR 916	High turbidity sample.	Filter visually turbid sample before analysis.
ERROR 918	Hardware or internal error.	The photometer needs service. Contact the distributor.

Symptom	Explanation	Action
HHH	1. Measured value exceeds 30.0 g/L (3.00 g/dL, 3,000 mg/dL, 1.86 mmol/L).	1a. Check expiry dates of the microcuvettes. 1b. For sample above 30.0 g/L (3.00 g/dL, 3,000 mg/dL, 1.86 mmol/L) use a suitable laboratory method.
No characters on the display	1. The photometer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order.	1a. Check that the power adapter is connected to the photometer and the AC power supply. 1b. Check that the cable is not damaged. 2. Turn off the photometer and replace the batteries, 5 alkaline type AA. 3. The photometer needs service. Contact the distributor.
The display contains erroneous characters	1. The display is out of order. 2. The microprocessor is out of order.	1. The photometer needs service. Contact the distributor. 2. The photometer needs service. Contact the distributor.
The display shows 'LOWBAT'	1. The batteries need to be replaced. 2. If on AC power, the power adapter or the circuit board is out of order.	1. Turn off the photometer and replace the batteries, 5 alkaline type AA. 2a. Check that the power adapter is properly connected and working. 2b. The photometer needs service. Contact the distributor.
The display does not switch from 'LHb' to 'READY' or from 'READY' to 'MEASURING'	1. The cuvette holder sensor is out of order.	1. The photometer needs service. Contact the distributor.

Symptom	Explanation	Action
Measurement on control materials is out of range – either too high or too low	1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The controls are beyond their expiry dates or have been improperly stored. 4. The control has not been mixed properly and/or is not at room temperature. 5. The microcuvette has not been placed in the photometer within 60 seconds of filling. 6. Air bubbles in the microcuvette. 7. The optical parts are dirty. 8. The control is not suitable for use with the HemoCue Plasma/Low Hb system. 9. The calibration of the photometer has been changed.	1. Check the expiry date and the storage conditions of the microcuvettes. 2. Remeasure the control with a new microcuvette. 3. Check the expiry date and the storage conditions of the control. Remeasure the control with a new microcuvette. If the problem continues, contact the manufacturer of the control. 4. Make sure that the control is mixed properly and at room temperature. If the problem continues, contact the manufacturer of the control. 5. Remeasure the control with a new microcuvette. 6. Check the microcuvette for air bubbles. Remeasure the control with a new microcuvette. 7. Clean the optical parts as described in the <i>Maintenance</i> section. 8. Contact the distributor for control information. 9. The photometer needs service. Contact the distributor.
Measurement on samples are higher or lower than anticipated	1. Improper sampling technique. 2. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 3. The optical eye of the microcuvette is contaminated. 4. The sample has not been mixed properly. 5. Air bubbles in the microcuvette. 6. The optical parts are dirty. 7. The calibration of the photometer has been changed.	1. See the <i>Measuring</i> section in this manual. 2. Check the expiry date and the storage conditions of the microcuvettes. 3. Remeasure the sample with a new microcuvette. 4. Make sure the sample is properly mixed. 5. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 6. Clean the optical parts as described in the <i>Maintenance</i> section. 7. The photometer needs service. Contact the distributor.

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Specifications

Intended Purpose

The HemoCue Plasma/Low Hb System is intended for quantitative determination of low levels of hemoglobin in plasma and serum specimens, aqueous solutions, or supernatant from erythrocyte suspensions using a specially designed photometer, the HemoCue Plasma/Low Hb Photometer and specially designed microcuvettes, the HemoCue Plasma/Low Hb Microcuvettes. The HemoCue Plasma/Low Hb System is an automated system intended for clinical laboratory testing by healthcare professionals and laboratory personnel. The HemoCue Plasma/Low Hb System is for *in vitro* diagnostic use only. The HemoCue Plasma/Low Hb Photometer is only to be used with HemoCue Plasma/Low Hb Microcuvettes.

Principle of the method

The reaction in the microcuvette is a modified azidemethemoglobin reaction. The erythrocytes are hemolyzed to release the hemoglobin. The hemoglobin is converted to methemoglobin and then combined with azide to form azidemethemoglobin. The measurement takes place in the photometer in which the transmittance is measured and the absorbance and hemoglobin level is calculated. The absorbance is directly proportional to the hemoglobin concentration.

Principle of the procedure

The system consists of a photometer together with the microcuvettes. The microcuvette serves both as a pipette and as a measuring cuvette and is for single-use only. A sample of approximately 20 µL is drawn into the cavity by capillary action. The photometer measures at two wavelengths in order to compensate for a certain degree of turbidity, and the hemoglobin level is calculated and presented (limitations apply as described in the paragraph *Limitations of the procedure*). The system is factory calibrated against hemoglobincyanide (HiCN) method, the international reference method for hemoglobin determination recommended by International Council for Standardization in Haematology (ICSH). No further calibration is needed.

Warning and precautions

The microcuvettes are for *in vitro* diagnostic use only. Always handle blood specimens with care as they may be infectious. Consult local environmental authorities for proper disposal.

Storage and handling

Operating temperature for the system is 15–30 °C (59–86 °F). The system should not be operated at < 5% or > 90% non-condensing humidity.

HemoCue Plasma/Low Hb Microcuvettes

The microcuvettes are to be stored at 15–30 °C (59–86 °F). Do not refrigerate. Use the microcuvettes prior to the expiration date that is printed on the package. Once the seal of the vial is broken, the microcuvettes are stable for three months. Keep the vial properly closed. All unused microcuvettes should remain in the original package.

HemoCue Plasma/Low Hb Photometer

The photometer can be stored and transported in temperatures between 0–50 °C (32–122 °F). Allow the photometer to reach operating temperature before use. The photometer should not be stored at < 5% or > 90% non-condensing humidity.

Specimen collection and preparation

Plasma/serum specimens and aqueous solutions in the form of irrigating fluid containing hemoglobin, may be used. Anticoagulant EDTA, Li-Heparin or Na-Heparin should be used for venous samples. Mix the specimen thoroughly before use. Supernatant from erythrocyte suspensions may be used. According to standardized laboratory procedures, centrifuge the erythrocyte suspension and carefully transfer the supernatant phase to a separate tube. Be cautious not to transfer any intact red cells from the precipitate. Do not let more than four hours pass between sampling and centrifugation. Visually turbid samples should be filtered before analysis using a 0.2 µm filter.

Materials required but not provided

- HemoCue Plasma/Low Hb Microcuvettes
- Protective gloves
- Pipette or other transfer device
- Lint-free wipe (non-fraying)
- Hydrophobic surface

Quality Control

If an external quality control test is required by local or other regulations, only use controls recommended by HemoCue, see relevant package insert for more information.

Measuring range

The measuring range is 0.3–30.0 g/L (0.03–3.00 g/dL, 30–3,000 mg/dL, 0.02–1.86 mmol/L).

Results

The measured hemoglobin value is read directly from the HemoCue Plasma/Low Hb Photometer. No calculations are necessary.

Limitations of the procedure

- If 'HHH' is displayed, the result exceeds the measuring range of the system.
- For values above 30.0 g/L (3.00 g/dL, 3,000 mg/dL, 1.86 mmol/L) use a suitable laboratory method.
- Sulfhemoglobin is not measured with this method.
- Levels of bilirubin (non-conjugated) up to 20 mg/dL and bilirubin (conjugated) up to 40 mg/dL do not influence the assay.
- Levels of intralipid up to 0.5 g/L do not influence the assay.
- The performance characteristics of this system have not been determined using samples obtained from uremic patients.

Specific performance characteristics

The results given below are from a HemoCue Plasma/Low Hb Photometer standardized against the HiCN reference method.

Linearity and Limit of Detection

The system is linear within displayed range 0–30.0 g/L (0–3.00 g/dL, 0–3,000 mg/dL, 0–1.86 mmol/L) (CLSI Document EP6-A).

Limit of detection (LoD) is 0.2 g/L (0.02 g/dL, 20 mg/dL, 0.01 mmol/L) (CLSI Document EP17-A2).

Repeatability and Within-laboratory precision

Results summarized below were determined for HemoCue Plasma/Low Hb System according to the CLSI Document EP5-A2.

	N	Mean (g/L)	Repeatability		Within-laboratory precision	
			SD (g/L)	CV %	SD (g/L)	CV %
Hemoglobin in 0.9% NaCl	400	1.4	0.05	3.4	0.05	3.6
	400	14.8	0.11	0.7	0.19	1.3
	400	24.7	0.17	0.7	0.28	1.1
Hemoglobin in plasma.	400	1.8	0.05	2.9	0.08	4.3
	400	14.8	0.14	1.0	0.21	1.4
	400	24.5	0.18	0.7	0.27	1.1

Correlation studies

Summarized below are the results of comparison studies between the HemoCue Plasma/Low Hb System and the HiCN reference method.

Study 1-2: Aqueous samples from surgical procedures.

Study 3-4: Erythrocyte suspensions.

Study	N	Min g/L	Max g/L	Regression line	Correlation coefficient (r)	Mean bias g/L
1	102	0.1	6.4	0.98 x HiCN – 0.05	0.98	–0.08
2	88	0.4	16.1	0.99 x HiCN + 0.01	0.99	–0.02
3	20	0.7	2.5	1.15 x HiCN – 0.02	0.99	0.19
4	64	0.4	6.5	1.02 x HiCN + 0.08	1.00	0.12

Technical Specifications

Dimensions: 160 x 210 x 90 mm (6.29 x 8.26 x 3.54 inches)

Weight (batteries included): 690g (1.5 pounds).

Battery: 5 batteries type AA, 1.5V Alkaline

Power adapter: CE marked

Only use the power adapter listed below:

Type	Input	Output
FW8002.1/12	100 V~240 V~/50–60 Hz/160–80 mA	12 V  /600 mA

Photometer electrical rating: 12 V , 0.5 W

Pollution degree: 2

Overvoltage category: II

Altitude: up to 2000 m above sea level

The instrument is made for continuous mode

IVD Medical Device Regulation

The HemoCue Plasma/Low Hb system complies with the IVD Medical Device Regulation (EU) 2017/746 and carries the CE mark.

EMC and Electrical Safety

The photometer has been tested and complies with following standards:

- IEC/EN 61010-1
- IEC/EN 61010-2-101
- IEC/EN 61326-2-6 (including applicable parts of IEC/EN 61326-1)

The photometer has been tested for indoor use.

It is recommended to evaluate the electromagnetic environment prior to use of the device.

Do not use this device in close proximity to sources of strong electromagnetic radiation (e.g unshielded intentional RF sources), as these can interfere with the proper operation.

Note: It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.

Note: It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.

Warranty

The photometer carries a 24-month warranty from the day of receipt. After the expiration date of the warranty, maintenance and repairs are offered at a fixed price. Any other use of the system than recommended by the manufacturer will void the warranty.

Service and disposal

The photometer should be cleaned as recommended under *Maintenance* prior to service or disposal. Dispose of the battery lid and cuvette holder as potentially infectious waste. Remove the batteries. Clean the photometer cover and optical parts according to section *Maintenance*. Disinfect the photometer external surfaces and optical parts according to local procedures. Wait for three weeks, then dispose of the photometer according to local procedures for electronic waste. For further information, contact HemoCue distributor or HemoCue AB.

Spare parts

Power adapter*
Battery lid
Cuvette holder

Optional items

HemoCue Cleaner

** Not available in all countries*

Reporting of serious incidents

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it to HemoCue AB or local distributor, and to your national authority as applicable according to local requirements.

Revision history

Rev. 230127 First revision to comply with IVDR

Symbols used

Caution

Consult instructions for use

In vitro diagnostic medical device

Biological risk

Manufacturer and date of manufacture

Serial number

Catalogue number

Temperature limitation

Humidity limitation

Only valid within the European Community. Indicates separate collection for waste of electrical and electronic equipment.

DC inlet

'OFF' (power)

'ON' (power)

Class II equipment
on power adapter

For indoor use only
on power adapter

References

- HemoCue Plasma/Low Hb package insert.
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A2
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline; CLSI Document EP6-A
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition; CLSI Document EP17-A2

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- Prospecto de HemoCue Plasma/Low Hb.
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A2
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline; CLSI Document EP6-A
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition; CLSI Document EP17-A2

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- HemoCue Plasma/Low Hb添付文書
- Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood; Approved Standard; CLSI Document H15
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline; CLSI Document EP5-A2
- Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline; CLSI Document EP6-A
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition; CLSI Document EP17-A2

© 2023 HemoCue AB

HemoCue AB | Kuvettgatan 1 | SE-262 71 Ängelholm | Sweden | Phone +46 77 570 02 10
info@hemocue.se | hemocue.com

900333 230922 EN ES JA